

MACHEREY-NAGEL

30 years of CHROMABOND®

Chromatography



Successful enrichment and purification

- Solid phase extraction (SPE) columns
- Established silica technology
- Innovative polymer-phases

 **Teknolab as**  
[www.teknolab.no](http://www.teknolab.no)

**MACHEREY-NAGEL**


[www.mn-net.com](http://www.mn-net.com)



# 30 years of CHROMABOND®

## Your advantages of 30 years of experience

- Robust standard silica phases
- Innovative RP- and mixed-mode phases
- Spherical polymer-based adsorbents



Find more than 4000 applications online at [www.mn-net.com/apps](http://www.mn-net.com/apps)

## Ordering information

Phase	Description	Volume	Adsorbent weight	Particle size	Pack of	REF
<b>The new one · N-vinylpyrrolidone-divinylbenzene resin</b>						
CHROMABOND® HLB	Hydrophilic-lipophilic balanced	3 mL	200 mg	60 µm	30	730924
CHROMABOND® HLB	Hydrophilic-lipophilic balanced	6 mL	200 mg	60 µm	30	730926
<b>The experts · Polystyrene-divinylbenzene resin</b>						
CHROMABOND® HR-X	Hydrophobic resin	3 mL	200 mg	85 µm	30	730931
CHROMABOND® HR-X	Hydrophobic resin	6 mL	200 mg	85 µm	30	730938
CHROMABOND® HR-XC	Strong mixed-mode cation exchanger	3 mL	200 mg	85 µm	30	730952
CHROMABOND® HR-XC	Strong mixed-mode cation exchanger	6 mL	150 mg	85 µm	30	730957
CHROMABOND® HR-XA	Strong mixed-mode anion exchanger	3 mL	200 mg	85 µm	30	730951
CHROMABOND® HR-XA	Strong mixed-mode anion exchanger	6 mL	150 mg	85 µm	30	730958
CHROMABOND® HR-XCW	Weak mixed-mode cation exchanger	3 mL	200 mg	85 µm	30	730739
CHROMABOND® HR-XCW	Weak mixed-mode cation exchanger	6 mL	150 mg	85 µm	30	730737
CHROMABOND® HR-XAW	Weak mixed-mode anion exchanger	3 mL	200 mg	85 µm	30	730748
CHROMABOND® HR-XAW	Weak mixed-mode anion exchanger	6 mL	150 mg	85 µm	30	730749
<b>The classics · Silica based phases</b>						
CHROMABOND® SiOH	Standard silica, polar	3 mL	500 mg	45 µm	50	730073
CHROMABOND® SiOH	Standard silica, polar	6 mL	1000 mg	45 µm	30	730075
CHROMABOND® C <sub>18</sub> ec	Modified encapped C <sub>18</sub> silica, nonpolar	3 mL	200 mg	45 µm	50	730012
CHROMABOND® C <sub>18</sub> ec	Modified encapped C <sub>18</sub> silica, nonpolar	6 mL	500 mg	45 µm	30	730014

Further filling amounts and hardware types (CHROMABOND® LV, CHROMAFIX®, CHROMABOND® MULTI 96) online at [www.mn-net.com](http://www.mn-net.com) or on request.

Vacuum manifold	REF
CHROMABOND Vacuum manifold, complete for 12 positions, consists of: glass cabinet with lid, lid gasket removable stainless steel needles on lower side of lid, vacuum gauge, control valve, valves and caps, PP-tank and variable rack	730150



Flyer Toolbox 30 years of CHROMABOND 2017 en1 / 10.2017



Company / distributor name  
Address  
Telephone  
E-Mail / website

[www.mn-net.com](http://www.mn-net.com)

# MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE / International:  
Tel.: +49 24 21 969-0  
Fax: +49 24 21 969-199  
E-mail: [info@mn-net.com](mailto:info@mn-net.com)

CH:  
Tel.: +41 62 388 55 00  
Fax: +41 62 388 55 05  
E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

FR:  
Tel.: +33 388 68 22 68  
Fax: +33 388 51 76 88  
E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

US:  
Tel.: +1 484 821 0984  
Fax: +1 484 821 1272  
E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

