

MACHEREY-NAGEL
CHROMAFIL[®] Xtra
Syringe filters

Chromatography



For reliable and successful filtrations

- PET syringe filters
- PTFE syringe filters
- PES syringe filters
- Polyamide / Nylon syringe filters

 **Teknolab as**
www.teknolab.no

MACHEREY-NAGEL
www.mn-net.com


Since 1911

CHROMAFIL® Xtra

CHROMAFIL® Xtra syringe filters

CHROMAFIL® syringe filters are used for the filtration of suspended matter from liquid samples (1–100 mL) or gases. With CHROMAFIL®, rapid purification and removal of particles is very simple: just place the filter on the syringe, and you are ready for filtration. Contamination of sensitive instrumentation by solid impurities can be avoided, thus increasing lifetime of chromatographic columns and equipment. The filters can be used for sample preparation in HPLC, GC, ICP, AAS, TOC, DOC, IR, NMR, photometry, spectroscopy, etc.

Product	Membrane pore size	Membrane dia.	Pack of	REF	Price
Polyester (PET), for polar as well as nonpolar solvents					
CHROMAFIL® Xtra PET-20/13	0.20 µm	13 mm	100	729222	
CHROMAFIL® Xtra PET-20/25	0.20 µm	25 mm	100	729221	
CHROMAFIL® Xtra PET-45/13	0.45 µm	13 mm	100	729223	
CHROMAFIL® Xtra PET-45/25	0.45 µm	25 mm	100	729220	
Polytetrafluoroethylene (PTFE), for nonpolar liquids and gases					
CHROMAFIL® Xtra PTFE-20/13	0.20 µm	13 mm	100	729208	
CHROMAFIL® Xtra PTFE-20/25	0.20 µm	25 mm	100	729207	
CHROMAFIL® Xtra PTFE-45/13	0.45 µm	13 mm	100	729209	
CHROMAFIL® Xtra PTFE-45/25	0.45 µm	25 mm	100	729205	
Polyethersulfone (PES), for aqueous and slightly organic liquids					
CHROMAFIL® Xtra PES-20/25	0.20 µm	25 mm	100	729240	
CHROMAFIL® Xtra PES-45/25	0.45 µm	25 mm	100	729241	
Polyamide / Nylon (PA), for aqueous organic and aqueous medium polar liquids					
CHROMAFIL® Xtra PA-20/13	0.20 µm	13 mm	100	729248	
CHROMAFIL® Xtra PA-20/25	0.20 µm	25 mm	100	729212	
CHROMAFIL® Xtra PA-45/13	0.45 µm	13 mm	100	729249	
CHROMAFIL® Xtra PA-45/25	0.45 µm	25 mm	100	729213	
Special filters for ion chromatography (IC), for aqueous liquids					
CHROMAFIL® Xtra IC-45/25	0.45 µm	25 mm	100	729258	



Company / distributor name
Address
Telephone
E-Mail / website



www.mn-net.com

MACHERY-NAGEL



MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

CH:
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

FR:
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

US:
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com

