

MACHEREY-NAGEL

CHROMABOND[®] QuEChERS

Chromatography



Convenient solutions for your pesticide residue analysis

- Proven CHROMABOND[®] Diamino (PSA)
- Pre-mixed in QuEChERS mixes



 **Teknolab as**
www.teknolab.no

MACHEREY-NAGEL

www.mn-net.com


Since 1911

CHROMABOND® QuEChERS

Ordering information

Extraction mixes

- 15 mL centrifuge tube with screw cap

Method	Mix No.	Content of mix	Volume	Quantity	REF	Price
EN 15662	Mix I	4 g MgSO ₄ , 1 g NaCl, 0.5 g Na ₂ H citrate · 1.5 H ₂ O, 1 g Na ₃ citrate · 2 H ₂ O	15 mL	50	730970	
AOAC 2007.01	Mix II	6 g MgSO ₄ , 1.5 g NaOAc	15 mL	50	730971	
Original	Mix XII	4 g MgSO ₄ , 1 g NaCl	15 mL	50	730648	

Clean-up mixes

- 15 mL centrifuge tube with screw cap
- 2 mL centrifuge tube with cap

Method	Mix No.	Content of mix	Volume	Quantity	REF	Price
EN 15662	Mix III	0.90 g MgSO ₄ , 0.15 g CHROMABOND® Diamino	15 mL	50	730972	
EN 15662	Mix IV	0.90 g MgSO ₄ , 0.15 g CHROMABOND® Diamino, 15 mg CHROMABOND® Carbon	15 mL	50	730973	
EN 15662	Mix V	0.90 g MgSO ₄ , 0.15 g CHROMABOND® Diamino, 45 mg CHROMABOND® Carbon	15 mL	50	730975	
EN 15662	Mix VI	0.90 g MgSO ₄ , 0.15 g CHROMABOND® Diamino, 150 mg CHROMABOND® C ₁₈ ec	15 mL	50	730974	
AOAC 2007.01	Mix XVII	0.15 g MgSO ₄ , 50 mg CHROMABOND® Diamino, 50 mg CHROMABOND® Carbon	2 mL	50	730996.2	
AOAC 2007.01	Mix XIX	0.15 g MgSO ₄ , 50 mg CHROMABOND® Diamino, 50 mg CHROMABOND® C ₁₈ ec	15 mL	50	730657	
AOAC 2007.01	Mix XX	1.20 g MgSO ₄ , 0.40 g CHROMABOND® Diamino	15 mL	50	730658	



Company / distributor name
Address
Telephone
E-Mail / website

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com

