

MACHEREY-NAGEL

Smart lab work...



Chromatography



...with vials and caps from MACHEREY-NAGEL

- High quality products for a reliable analysis
- Smart savings for selected products
- Vials and caps for HPLC and IC



MACHEREY-NAGEL

www.mn-net.com



Products for HPLC / Ion chromatography

Glass and plastic vials N 9 and N 11, inserts and closures

	Screw neck vial N 9 702283 1.5 mL 11.6 x 32 mm clear, label + scale	Screw neck vial N 9 702284 1.5 mL 11.6 x 32 mm amber, label + scale	Screw closure N 9 702027 N 9 PP bonded screw cap, blue, Silicone beige / PTFE white, slit	Screw closure N 9 702288.1 N 9 PP screw cap, blue, Silicone white / PTFE blue, slit	Screw closure N 9 702287.1 N 9 PP screw cap, blue, Silicone white / PTFE red
	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.
	Snap ring vial N 11 702713 1.5 mL 11.6 x 32 mm clear, label + scale	Snap ring vial N 11 702712 1.5 mL 11.6 x 32 mm amber, label + scale	Snap ring closure N 11 702064 N 11 PE snap ring cap, transparent, Silicone white / PTFE blue, cross-slit	Snap ring closure N 11 702710 N 11 PE snap ring cap, transparent, Silicone white / PTFE red	Snap ring closure N 11 702731 N 11 PE snap ring cap, transparent, Red Rubber / FEP colorless
	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.
	Insert for a.m. vials 702818 0.1 mL 5.7 x 29 mm with plastic spring	Insert for a.m. vials 702813 0.2 mL 6 x 31 mm 15 mm tip	Insert for a.m. vials 702825 0.3 mL 6 x 31 mm flat bottom	Micro-vial N 9 702135 PP screw neck vial with integrated 0.2 mL glass insert	Micro-vial N 11 702134 PP snap ring vial with integrated 0.2 mL glass insert
	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.
	PP screw neck vial N 9 702500 1.5 mL 11.6 x 32 mm transparent, filling lines	PP screw neck vial N 9 702009 0.3 mL 11.6 x 32 mm transparent, inner cone	PP screw neck vial N 9 702010 0.7 mL 11.6 x 32 mm transparent, round bottom insert	PP snap ring vial N 11 702809 0.3 mL 11.6 x 32 mm transparent, inner cone	PP snap ring vial N 11 702174 0.7 mL 11.6 x 32 mm transparent, round bottom insert
	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.	€ XX.XX per pk.

All products in packaging units of 100 pcs.

Accessories

Vial rack

702502
 50 position PP vial rack blue for all vials
 11.6 x 32 mm with flat bottom
 Dimensions: 190 x 100 x 22 mm, stackable

€ XX.XX per pc.

Vial container

702514
 81 position container blue for screw neck vials N 9 and snap ring vials N 11,
 130 x 130 x 45 mm, coded,
 with transparent lid (suitable for freezers)

€ XX.XX per pc.



Company / distributor name
 Address
 Telephone
 E-Mail / website

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com

