

# MACHEREY-NAGEL

## CHROMABOND<sup>®</sup> HLB

Chromatography

Request  
your free  
sample



Enhanced retention for polar compounds

- Hydrophilic-lipophilic balanced phase
- Wide range of applications – biological, pharmaceutical or environmental samples

Logo

**MACHEREY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



# CHROMABOND® HLB

## Advantages of CHROMABOND® HLB

- Applicable for a wide range of analyte polarities
- Enhanced retention for polar compounds
- High loadability and outstanding performance
- Water wettable – even if bed runs dry, SPE can be continued
- The alternative to Oasis® HLB

## Technical data

Hydrophilic-lipophilic balanced *N*-vinylpyrrolidone-divinylbenzene copolymer


Particle shape: spherical

pH stability: 1–14

Particle size: 60 µm and 30 µm

Pore size: 65 Å

Specific surface: 750 m<sup>2</sup>/g



## Ordering information

REF	Product	Particle size	Volume	Adsorbent weight	Pack of	Price	Unit
730921	CHROMABOND® HLB	60 µm	1 mL	30 mg	30		EUR
730921P30	CHROMABOND® HLB	30 µm	1 mL	30 mg	30		EUR
730922	CHROMABOND® HLB	60 µm	1 mL	100 mg	30		EUR
730922P30	CHROMABOND® HLB	30 µm	1 mL	100 mg	30		EUR
730923	CHROMABOND® HLB	60 µm	3 mL	60 mg	30		EUR
730923P30	CHROMABOND® HLB	30 µm	3 mL	60 mg	30		EUR
730924	CHROMABOND® HLB	60 µm	3 mL	200 mg	30		EUR
730924P30	CHROMABOND® HLB	30 µm	3 mL	200 mg	30		EUR
730944	CHROMABOND® HLB	60 µm	6 mL	150 mg	30		EUR
730944P30	CHROMABOND® HLB	30 µm	6 mL	150 mg	30		EUR
730926	CHROMABOND® HLB	60 µm	6 mL	200 mg	30		EUR

Further filling amounts and hardware types (CHROMABOND® LV, CHROMAFIX®, CHROMABOND® MULTI 96) online at [www.mn-net.com](http://www.mn-net.com) or on request.

Registered trademarks:


Oasis® Waters Corp. (USA)

CHROMABOND® MACHEREY-NAGEL GmbH & Co. KG (Germany)

CHROMAFIX® MACHEREY-NAGEL GmbH & Co. KG (Germany)



Teknolab AS  
Verkstedveien 29, 1424 Ski  
tlf. 66813470  
mail@teknolab.no / www.teknolab.no



Visit our application database at  
[www.mn-net.com/apps](http://www.mn-net.com/apps)

[www.mn-net.com](http://www.mn-net.com)

# MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: [info@mn-net.com](mailto:info@mn-net.com)

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-Mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

